



Online store [WWW.LYNXGEAR.LV](http://WWW.LYNXGEAR.LV)  
Company: SIA LATLYNX  
info@lynxgear.lv, +371 26159668  
Reg. No. 45403038162, VAT Reg. No.  
LV45403038162,  
Legal address: Baznicas iela 8,  
Aizkraukle, Aizkraukles novads,  
LV-5101, LATVIJA  
Bank Account: LV87UNLA0050021387494  
AS SEB Bank, SWIFT: UNLALV2X

---

## ORDER CONFIRMATION

---

Not Rated Yet

**Sales price 21,95 €**

Salesprice with discount

Incl. VAT 21%: 3,81 €



<b>Needle Gauge</b>	14G
<b>Total Length</b>	14.3 cm
<b>Needle Length</b>	8.8 cm
<b>Needle Diameter</b>	1.15 mm
<b>Weight (with shell)</b>	14 g
<b>Material</b>	Stainless steel
<b>Shelf Life</b>	5 years from manufacturing date
<b>Product Code</b>	H5000-4012D
<b>UDI</b>	06974582640289
<b>NSN</b>	6515-01-541-0635
<b>Certification</b>	CE-marked medical device

### Typical Uses

- Tension pneumothorax decompression
- Tactical and military medical kits
- Advanced life support (ALS) kits
- Prehospital emergency care

### Important Information

This product is intended for use by **trained medical professionals only.**



Online store [WWW.LYNXGEAR.LV](http://WWW.LYNXGEAR.LV)  
Company: SIA LATLYNX  
info@lynxgear.lv, +371 26159668  
Reg. No. 45403038162, VAT Reg. No.  
LV45403038162,  
Legal address: Baznicas iela 8,  
Aizkraukle, Aizkraukles novads,  
LV-5101, LATVIJA  
Bank Account: LV87UNLA0050021387494  
AS SEB Bank, SWIFT: UNLALV2X

---

## ORDER CONFIRMATION

---

Incorrect use may result in serious injury or death.

### About the Manufacturer

**CVN Medical Solutions** develops advanced prehospital care products focused on bleeding control, hypothermia prevention and tactical medical operations. Their equipment supports both civilian emergency services and military applications.

Beyond medical devices, CVN Medical Solutions also provides certified **TCCC** and **TECC** training programs. The company operates under strict quality standards including **ISO 9001**, **ISO 14001**, **ISO 13485** and **AQAP**, ensuring high reliability and compliance.

### FAQ

**What is a tension pneumothorax needle used for?** For air trapped in the chest cavity (tension pneumothorax).

**Who can use this product?** For trained medical personnel only.

**Is this product sterile?** Sterile and ready for immediate use.

**Is this included in standard first aid kits?** Tactical trauma kits.

### Glossary

- **Tension pneumothorax:** A life-threatening condition caused by trapped air in the chest cavity.
- **Decompression:** The release of pressure from the chest to restore breathing.
- **Gauge (G):** Measurement of needle diameter.

### Safety Information

Improper use can result in severe injury or fatal complications. Use only in accordance with medical training and protocols.

### Manufacturer & GPSR Contact

**CVN Medical Oy**

Yrittäjätie 2

FI-60100 Seinäjoki

Finland

Email: This email address is being protected from spambots. You need JavaScript enabled to view it.

```
document.getElementById('cloakfd4e0ef72cc3404eee2de21925a78d0d').innerHTML = ''; var prefix = 'ma' + 'il' + 'to'; var path = 'hr' + 'ef' + '='; var addyfd4e0ef72cc3404eee2de21925a78d0d = 'info' + '@'; addyfd4e0ef72cc3404eee2de21925a78d0d = addyfd4e0ef72cc3404eee2de21925a78d0d + 'cvn' + '.' + 'fi'; var addy_textfd4e0ef72cc3404eee2de21925a78d0d = 'info' + '@' + 'cvn' + '.' + 'fi'; document.getElementById('cloakfd4e0ef72cc3404eee2de21925a78d0d').innerHTML += ''+addy\_textfd4e0ef72cc3404eee2de21925a78d0d+';
```