



Online store WWW.LYNXGEAR.LV
Company: SIA LATLYNX
info@lynxgear.lv, +371 26159668
Reg. No. 45403038162, VAT Reg. No.
LV45403038162,
Legal address: Baznicas iela 8,
Aizkraukle, Aizkraukles novads,
LV-5101, LATVIJA
Bank Account: LV87UNLA0050021387494
AS SEB Bank, SWIFT: UNLALV2X

ORDER CONFIRMATION

Not Rated Yet

Sales price 21,95 €

Salesprice with discount

Incl. VAT 21%: 3,81 €



Needle Gauge	14G
Total Length	14.3 cm
Needle Length	8.8 cm
Needle Diameter	1.15 mm
Weight (with shell)	14 g
Material	Stainless steel
Shelf Life	5 years from manufacturing date
Product Code	H5000-4012D
UDI	06974582640289
NSN	6515-01-541-0635
Certification	CE-marked medical device

Typical Uses

- Tension pneumothorax decompression
- Tactical and military medical kits
- Advanced life support (ALS) kits
- Prehospital emergency care

Important Information

This product is intended for use by **trained medical professionals only.**



Online store WWW.LYNXGEAR.LV
Company: SIA LATLYNX
info@lynxgear.lv, +371 26159668
Reg. No. 45403038162, VAT Reg. No.
LV45403038162,
Legal address: Baznicas iela 8,
Aizkraukle, Aizkraukles novads,
LV-5101, LATVIJA
Bank Account: LV87UNLA0050021387494
AS SEB Bank, SWIFT: UNLALV2X

ORDER CONFIRMATION

Incorrect use may result in serious injury or death.

About the Manufacturer

CVN Medical Solutions develops advanced prehospital care products focused on bleeding control, hypothermia prevention and tactical medical operations. Their equipment supports both civilian emergency services and military applications.

Beyond medical devices, CVN Medical Solutions also provides certified **TCCC** and **TECC** training programs. The company operates under strict quality standards including **ISO 9001**, **ISO 14001**, **ISO 13485** and **AQAP**, ensuring high reliability and compliance.

FAQ

What is a pneumothorax needle used for? by air trapped in the chest cavity (tension pneumothorax).

Who can use this product? for trained medical personnel only.

Is this product sterile? ile and ready for immediate use.

Is this included in standard first aid kits? tactical trauma kits.

Glossary

- **Tension pneumothorax:** A life-threatening condition caused by trapped air in the chest cavity.
- **Decompression:** The release of pressure from the chest to restore breathing.
- **Gauge (G):** Measurement of needle diameter.

Safety Information

Improper use can result in severe injury or fatal complications. Use only in accordance with medical training and protocols.

Manufacturer & GPSR Contact

CVN Medical Oy

Yrittäjäntie 2

FI-60100 Seinäjoki

Finland

Email: This email address is being protected from spambots. You need JavaScript enabled to view it.

```
document.getElementById('cloakda2d9e3f41fa8b7a50cc4b60c59a250e').innerHTML = ''; var prefix = 'ma' + 'il' + 'to'; var path = 'hr' + 'ef' + '='; var addyda2d9e3f41fa8b7a50cc4b60c59a250e = 'info' + '@'; addyda2d9e3f41fa8b7a50cc4b60c59a250e = addyda2d9e3f41fa8b7a50cc4b60c59a250e + 'cvn' + '.' + 'fi'; var addy_textda2d9e3f41fa8b7a50cc4b60c59a250e = 'info' + '@' + 'cvn' + '.' + 'fi';document.getElementById('cloakda2d9e3f41fa8b7a50cc4b60c59a250e').innerHTML += ''+addy\_textda2d9e3f41fa8b7a50cc4b60c59a250e+';
```